

TP N°03

Acid-base dosage

Introduction

➤ One of the first definitions of an acid and a base was proposed in 1887 by Arrhenius and Ostwald. According to them, an acid was a compound with mobile hydrogen which released H^+ protons into water while a base was a compound which released OH^- hydroxide ions in an aqueous medium.


➤ In 1923, Bronsted and Lowry questioned the definition of Arrhenius and Ostwald and proposed a new definition. They defined an acid as a substance capable of giving up a proton and a base as a substance capable of capturing a proton.

➤ At the same time, during the same year, Lewis proposed another definition. He defined an acid as an acceptor of electron pairs and a base as a donor of electron pairs.

➤ An acid-base couple consists of an acid and its conjugate base (Acid/Base).

➤ An acid-base reaction involves two acid-base couples: the Acid1/Base1 couple and the Acid2/Base2 couple.

➤ An acid-base reaction is a chemical transformation between the acid of one pair and the base of another acid/base pair, via an exchange of H^+ ions. The full equation is a linear combination of the two specific half-equations of each pair.



Principle of a dosage

A solution contains a dissolved chemical species A. **Determining this chemical species** means determining its **quantity of matter** or its **C_A concentration** in the solution. To measure A, A is reacted with a body B contained in a **solution of known concentration C_B** . The dosage reaction must be rapid, complete, easily observable.

- An acid-base dosage can be followed by:
 - **pH-metry**: we follow the evolution of the pH during the reaction.
 - **Colorimetry**: we use a colored indicator

A colored indicator is a reagent whose color depends on the medium (or pH).

It can be used to mark the end of a dosage if equivalence is reached in its turning zone.

Examples of colored indicators:

Indicator	Acid / tint	pH range (turning zone)	Base / tint
Helianthin	Red	3,1 -4,4	Yellow
Methyl red	Red	4,08 -6,0	Yellow
Bromothymole blue	Yellow	6,0 -7,6	Blue
Phenol-phthalein		8,2 -10,0	
Alizarin Yellow	Yellow	10,1 -12,2	Red

Objective:

Determine the molar concentration of an acid solution, using the colorimetric dosage.

Material:


Graduated or volumetric pipettes (10 mL), Suction device, Burette (25 or 50 mL), Erlenmeyer flask (100 mL), beaker (x2).

Operating mode:

1. Dosage of a strong acid with a strong base

The dosage of the **hydrochloric acid** will be carried out using a **sodium hydroxide** solution (NaOH) with a molar concentration $C_B = 0.1 \text{ mol/L}$ in the presence of phenolphthalein.

The reaction equation is:



Rapid dosing (determination of a framework for the equivalent volume V_e)

- Check that the burette stopcock is closed.
- Rinse the graduated burette with the titrant solution (NaOH) of a precise molar concentration ($C_B = 0.1 \text{ mol/L}$), then fill it.
- Adjust the liquid level to the zero level of the burette by draining the excess sodium hydroxide solution into the labeled beaker.
- Pour approximately 40 mL of solution S1 into a labeled beaker.
- Introduce into a 100 mL Erlenmeyer flask:
 - 10 mL of hydrochloric acid solution taken using a clean volumetric pipette fitted with a suction device,
 - 1 to 3 drops of phenolphthalein,
- Place the Erlenmeyer flask under the burette, shake manually without using the magnetic stirrer
- Add the titrant solution (mL per mL) and note the color of the solution by completing table 1:




Table 1:

V sol. titrant (ml)	1	2	3	4	5	6	7	8
Color								

Req.: The solution changed color when you added the equivalent volume of titrant solution (V_e).

Indicate approximately this volume (by a frame); $V_1 \text{ mL} < V_e < V_2 \text{ mL}$

Precise dosage (known as “drop dosage”)

Add the titrant solution until the color changes (equivalence point) while respecting the following instructions:

- quickly at the beginning, a Volume $\leq V_1$
- then drop by drop as the color change approaches (equivalence point).
- Read the equivalent volume and note the color of the solution by completing table 2.
- Repeat the operation two to three times

Table2:

	1st test	2nd test	3rd test
V_e (ml)			
Color			

- Deduce the average equivalence volume (Avg. V_{eq})?
- Calculate the concentration of the hydrochloric acid solution to be titrated S2?
- Deduce the concentration of the solution S1

2. Determination of a diacid using a strong base

The dosage of sulfuric acid (H_2SO_4) will be carried out using a sodium hydroxide solution with a molar concentration $C_B = 0.1 \text{ mol/L}$ in the presence of phenolphthalein.

Rapid dosage (determination of a framework for the equivalent volume V_e)

- Check that the burette stopcock is closed.
- Rinse the graduated burette with the titrant solution (KOH or NaOH) of very precise molar concentration ($C_B = 0.1 \text{ mol/L}$), then fill it.
- Adjust the liquid level to the zero level of the burette by draining the excess soda (potash) solution into the beaker labeled (Recovery of used products).
- Pour approximately 50 mL of the sulfuric acid solution into a labeled beaker.
- Introduce into a 100 mL Erlenmeyer flask:
 - 10 mL of sulfuric acid solution taken using a clean volumetric pipette fitted with a propipette,
 - 3 drops of phenolphthalein,
- Place the Erlenmeyer flask under the burette. Stir the solution manually without using the magnetic stirrer.
- Add the titrant solution (mL per mL) and note the color of the solution by completing

table 3.

Table 3:

V sol. titrant (ml)	1	2	3	4	5	6	7	8
Color								

Req.: The solution changed color when you added the equivalent volume of titrant solution (V_e).

Indicate approximately this volume (by a frame); **$V_1 \text{ mL} < V_e < V_2 \text{ mL}$**

Precise dosage (known as “drop dosage”)

Add the titrant solution until the color changes (equivalence point) while respecting the following instructions:

- quickly at the beginning, a Volume $\leq V_1$
- then drop by drop as the color change approaches (equivalence point).
- Read the equivalent volume and note the color of the solution by completing table 4.
- Repeat the operation two to three times

Table 4:

	1st test	2nd test	3rd test
V_e (ml)			
Color			

- Deduce the average equivalence volume (Avg. V_{eq})?
- Calculate the concentration of the sulfuric acid solution to be titrated?
- What is the difference between first dosage and second dosage?

V. Questions:

1. Are there any spectator chemical species in both dosages? Which ones?
2. What are the acid/base pairs involved in each dosage?
3. Write the associated proton half-equations in the two dosages?
4. Derive the equation for the reaction of the dosage in both cases?
5. Give the definition of dosage equivalence? Deduce a relationship between the quantities of matter of oxonium ions (H_3O^+) (n_A) and hydroxide ions (OH^-) (n_B)?
6. List the chemical species present in the dosing beaker:
 - for a volume poured less than V_e^B ?
 - for a poured volume equal to V_e^B . What should the pH be at equivalence?
 - for a poured volume greater than V_e^B ?